

TUCAM-API

Guide for Development

Fuzhou Xintu Photonics Co., Ltd.

Copyright (c) 2011-2016 Xintu Photonics Co., Ltd.

(TUCESN)

All rights reserved.

Table of Contents

Table of Contents	1
1. Read before use.....	4
2. Brief Introduction	5
3. Overview	6
3.1 Layer Structure	6
3.2 Principle.....	6
3.3 Interface Type	6
3.4 Terminology.....	7
3.4.1 Capture Mode:	7
3.4.2 Picture Cell:	7
3.4.3 Trigger Mode:	7
3.4.4 Camera Status:	8
3.5 Interface List	8
4. Application Calls TUCAM-API.....	10
4.1 Initialization and Termination of Procedures	10
4.1.1 Interface:	10
4.1.2 Calling Sequence:	10
4.1.3 Driver Initialization:	10
4.1.4 Camera Initialization:.....	11
4.1.5 Camera Product Information:	11
4.1.6 Termination Procedure:	12
4.1.7 Sample Code:	12
4.2 Performance Obtaining and Setting	13
4.2.1 Interface:	13
4.2.2 Calling Sequence:	13
4.2.3 Performance Index:	13
4.2.3 Sample Code:	13
4.3 Property Extraction and Settings.....	15
4.3.1 Interface:	15
4.3.2 Calling Sequence:	15
4.3.3 Property Index:	15
4.3.3 Sample Code:	16
4.4 Memory Management	17
4.4.1 Interface:	17
4.4.2 Calling Sequence:	17
4.4.3 Frame Structure:	17
4.4.4 Sample Code:	18
4.5 Capture Control	20
4.5.1 Interface:	20
4.5.2 Calling Sequence:	21
4.5.2 Capture Mode Index:.....	21
4.5.3 Sample Code:	21
4.6 File Control	22

4.6.1 Interface:	22
4.6.2 Calling Sequence:	22
4.6.3 File Structure:	22
4.6.4 Sample Code:	23
4.7 Extending Control	24
4.7.1 Interface:	24
4.7.2 Calling Sequence:	24
4.7.3 File Structure:	25
4.7.4 Sample Code:	25
5. Reference	27
5.1 Types and Constants	27
5.1.1 TUCAMRET Error code:	27
5.1.2 TUCAM_IDINFO Product Information Code:	28
5.1.3 TUCAM_IDCAPA Performance Code:	29
5.1.4 TUCAM_IDPROP Property Code:	29
5.1.5 TUCAM_CAPTURE_MODES Capture Mode Code:	29
5.1.6 TUIMG_FORMATS Image Format Code:	30
5.1.7 TUREG_TYPE Register Type Code:	30
5.1.8 TUCAM_TRIGGER_EXP Trigger Exposure Mode Code:	30
5.1.9 TUCAM_TRIGGER_EDGE Trigger Excitation Edge Code:	30
5.1.10 TUFRM_FORMATS Frame Format Code:	30
5.2 Structure	31
5.2.1 Initialization:	31
5.2.2 Open Camera:	31
5.2.3 Camera Information:	31
5.2.4 Performance/Property Value Text:	31
5.2.5 Performance Property:	32
5.2.6 Property Attribute:	32
5.2.7 ROI Attribute:	32
5.2.8 Trigger Attribute:	32
5.2.9 Frame Structure:	33
5.2.10 File Saving:	33
5.2.11 Video Recording Saving:	34
5.2.12 Read/Write of Register	34
5.3 Functions	35
TUCAM_Api_Init	35
TUCAM_Api_Uninit	35
TUCAM_Dev_Open	36
TUCAM_Dev_Close	37
TUCAM_Dev_GetInfo	37
TUCAM_Capa_GetAttr	38
TUCAM_Capa_GetValue	39
TUCAM_Capa_SetValue	40
TUCAM_Capa_GetValueText	40

TUCAM_Prop_GetAttr	41
TUCAM_Prop_GetValue	42
TUCAM_Prop_SetValue	43
TUCAM_Prop_GetValueText	44
TUCAM_Buf_Alloc	44
TUCAM_Buf_Release	45
TUCAM_Buf_AbortWait	46
TUCAM_Buf_WaitForFrame	47
TUCAM_Buf_CopyFrame	48
TUCAM_Cap_SetROI	49
TUCAM_Cap_GetROI	50
TUCAM_Cap_SetTrigger	50
TUCAM_Cap_GetTrigger	53
TUCAM_Cap_DoSoftwareTrigger	53
TUCAM_Cap_Start	54
TUCAM_Cap_Stop	55
TUCAM_File_SaveImage	56
TUCAM_Rec_Start	56
TUCAM_Rec_AppendFrame	57
TUCAM_Rec_Stop	58
TUCAM_Reg_Read	59
TUCAM_Reg_Write	59

1. Read before use

The document and software example code are the internal document and release content of TUCSEN, which are provided for the purpose of enabling the user to create and use the Apps of TUCSEN Digital Camera.

This document and software sample code are disclosed for the above purpose only, which shall not constitute the owner's permission, transfer or any other power.

All risks and the results of using software documentation still depend on the user.

This document may contain technical or typographical errors, and it does not guarantee any damage resulting from such mistakes or text.

TUCSEN does not make commitment to update or to keep the information contained in the current document.

All brands and product names are the trademarks or registered trademarks of their respective owners.

TUCSEN reserves all rights of the document copyright.

In the absence of prior written permission of TUCSEN, any part of the document can not be reproduced, transmitted, transcribed or stored in a retrieval system or translated into any language or computer language, in any form or by any means and any methods such as electronic, mechanical, electromagnetic, optical, chemical, manual or other ones.

2. Brief Introduction

This manual describes the detailed specification of TUCAM-API as how to use TUCSEN Digital Camera. TUCAM-API Software Development Kit is called "SDK". The part of TUCAM-API controls the digital camera is called "module".

SDK contains source code modules and one sample application that demonstrates how to access TUCAM-API. SDK user can freely use the software in any way they like, such as partial modification of the source code or the creation of a completely separate project.


It is particularly easy to understand the SDK design. For this reason, it limits the number of function interfaces to a minimum, and the function call format is written in C language.

Parts of the extended functions are the additional features which can be used for certain digital cameras.

The values of the digital cameras may be different, which depend on the models of digital cameras used to capture images. The values should be simply seen as the guideline rather than the exact values.

3. Overview

3.1 Layer Structure



TUCSEN Digital Camera is connected with the drivers of digital camera of the different operating systems via SDK to achieve the functions of controlling the digital camera and capture of images and data.

The SDK currently supports only Windows systems.

3.2 Principle

The specific bus interfaces and libraries of digital cameras are packaged via TUCAM-API. You only need to access the TUCAM-API Layer. The Module Layer provides more advanced TUCAM-API integration. Modules can be constantly updated to access to the new cameras and provide new interface technology without recompiling of your software.

3.3 Interface Type

TUCAM-API functions can be divided into many types:

- Start/End Processing;
- Camera information collection;
- Performance/property extraction and settings;
- Memory management;
- Capture control;
- Document control;

Extending control.

TUCAM-API does not include the program for displaying images. Since it is difficult to predict some methods of displaying images, which depends on the application and it is impossible to support all of these common modules. When calling the display program, it will detect whether the image is updated, and it will draw the image after the update. For more detailed information, please refer to the sample source code.

3.4 Terminology

3.4.1 Capture Mode:

Camera capture modes are divided into the following 2 categories:

Sequence mode (stream mode): it is used to capture continuous image data.

Trigger mode: The camera captures the images by the external signals. We call this option as "trigger mode", and you can call `TUCAM_Cap_SetTrigger()` to configure this option. We also call the external signal as "external trigger".

3.4.2 Picture Cell:

It is generally two-dimensional, with vertical and horizontal directions.

Frame: it is a unit of image data. For one frame, the pixel data is aligned from left to right and from top to bottom. This is a series of image data unit.

3.4.3 Trigger Mode:

Standard Mode: When the camera receives the level signal (determined by the activation edge), it will start the image capturing of one frame or multiple frames. The number of captured frames is determined by the configuration parameters. Refer to `TUCAM_TRIGGER_ATTR` Structure.

Synchronization Mode: When the camera receives the level signal (determined by the activation edge), it will begin the exposure; after receiving the opposite level signal, it will end the exposure, and start the capture of image data. That is, the completion of exposure and read-out of each frame is completely synchronized with the external trigger signal.

Global Mode: Pre-trigger should be conducted prior to triggering the camera. When the camera receives a level signal (determined by the activation edge) or the exposure time set for the software, the reset operation currently in progress will be ended. When the exposure is ended, the image data will be captured, and the pre-trigger will be restarted. This mode is used

to control the camera of shutter exposure mode to achieve global exposure mode.

Exposure Mode:

Exposure Time: after receiving the trigger signal, it is determined by the exposure time set by TUIDP_EXPOSURETM.

Level Width: after receiving the trigger signal, the exposure time is determined by the level width.

NOTE: The Standard Mode and Global Mode can configure both options while Synchronization Mode can only be the level width.

Types of Excitation Level:

Rising Edge: The trigger level starts exposure at the rising edge

Falling Edge: The trigger level starts exposure at the falling edge

3.4.4 Camera Status:

The camera status determines which functions can be called. Some functions will change the state of the camera. Four kinds of camera statuses are described as follows:

Unstable: parameter setting and call of other functions, but they are not in the status to be set.

Stable: the parameters and functions are set, but the image capture cannot be started as there is no frame memory created.

Preparation: frame memory has been created, and the image capture can be started.

Busy: image capturing is being executed.

3.5 Interface List

// Initialize uninitialized and misc.

TUCAMRET TUCAM_Api_Init(PTUCAM_INIT pInitParam);

TUCAMRET TUCAM_Api_Uninit ();

TUCAMRET TUCAM_Dev_Open (PTUCAM_OPEN pOpenParam);

TUCAMRET TUCAM_Dev_Close (HDTUCAM hTUCam);

// Get some device information (VID/PID/Version)

TUCAMRET TUCAM_Dev_GetInfo (HDTUCAM hTUCam, PTUCAM_VALUE_INFO pInfo);

// Capability control see enumerate TUCAM_IDCAPA

TUCAMRET TUCAM_Capa_GetAttr (HDTUCAM hTUCam, PTUCAM_CAPA_ATTR pAttr);

TUCAMRET TUCAM_Capa_GetValue (HDTUCAM hTUCam, INT32 nCapa, INT32 *pnVal);

```

TUCAMRET  TUCAM_Capa_SetValue (HDTUCAM hTUCam, INT32 nCapa, INT32 nVal);
TUCAMRET  TUCAM_Capa_GetValueText (HDTUCAM hTUCam, PTUCAM_VALUE_TEXT pVal);

// Property control see enumerate PTUCAM_PROP_ATTR
TUCAMRET  TUCAM_Prop_GetAttr (HDTUCAM hTUCam, PTUCAM_PROP_ATTR pAttr);
TUCAMRET  TUCAM_Prop_GetValue (HDTUCAM hTUCam, INT32 nProp, DOUBLE *pdbVal, INT32 nChn);
TUCAMRET  TUCAM_Prop_SetValue (HDTUCAM hTUCam, INT32 nProp, DOUBLE dbVal, INT32 nChn );
TUCAMRET  TUCAM_Prop_GetValueText (HDTUCAM hTUCam, PTUCAM_VALUE_TEXT pVal, INT32
nChn );

// Buffer control
TUCAMRET  TUCAM_Buf_Alloc (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
TUCAMRET  TUCAM_Buf_Release (HDTUCAM hTUCam);
TUCAMRET  TUCAM_Buf_AbortWait (HDTUCAM hTUCam);
TUCAMRET  TUCAM_Buf_WaitForFrame (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
TUCAMRET  TUCAM_Buf_CopyFrame (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);

// Capturing control
// ROI
TUCAMRET  TUCAM_Cap_SetROI (HDTUCAM hTUCam, TUCAM_ROI_ATTR roiAttr);
TUCAMRET  TUCAM_Cap_GetROI (HDTUCAM hTUCam, PTUCAM_ROI_ATTR pRoiAttr);
// Trigger
TUCAMRET  TUCAM_Cap_SetTrigger (HDTUCAM hTUCam, TUCAM_TRIGGER_ATTR tgrAttr);
TUCAMRET  TUCAM_Cap_GetTrigger (HDTUCAM hTUCam, PTUCAM_TRIGGER_ATTR pTgrAttr);
TUCAMRET  TUCAM_Cap_DoSoftwareTrigger(HDTUCAM hTUCam);    // in trigger mode

// Capturing
// uiMode see enumerate TUCAM_CAPTURE_MODES
TUCAMRET  TUCAM_Cap_Start(HDTUCAM hTUCam, UINT32 uiMode);
TUCAMRET  TUCAM_Cap_Stop (HDTUCAM hTUCam);

// File control
// Image
TUCAMRET  TUCAM_File_SaveImage (HDTUCAM hTUCam, TUCAM_FILE_SAVE fileSave);
// Video
TUCAMRET  TUCAM_Rec_Start(HDTUCAM hTUCam, TUCAM_REC_SAVE recSave);
TUCAMRET  TUCAM_Rec_AppendFrame(HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
TUCAMRET  TUCAM_Rec_Stop (HDTUCAM hTUCam);

// Extended control
TUCAMRET  TUCAM_Reg_Read (HDTUCAM hTUCam, TUCAM_REG_RW regRW);
TUCAMRET  TUCAM_Reg_Write(HDTUCAM hTUCam, TUCAM_REG_RW regRW);

```

4. Application calls TUCAM-API

When TUCAM-API is being used to control camera, the function call should be executed in accordance with the following calling procedure:

- ◆ Initialize camera
- ◆ Set camera parameters
- ◆ Start capturing data
- ◆ Make sure that shooting has been completed and the data have been obtained
- ◆ Perform camera termination processing

4.1 Initialization and Termination of Procedures

4.1.1 Interface:

// Initialize uninitialized and misc.

```
TUCAMRET TUCAM_Api_Init(PTUCAM_INIT pInitParam);
```

```
TUCAMRET TUCAM_Api_Uninit ();
```

```
TUCAMRET TUCAM_Dev_Open (PTUCAM_OPEN pOpenParam);
```

```
TUCAMRET TUCAM_Dev_Close (HDTUCAM hTUCam);
```

// Get some device information (VID/PID/Version)

```
TUCAMRET TUCAM_Dev_GetInfo (HDTUCAM hTUCam, PTUCAM_VALUE_INFO pInfo);
```

4.1.2 Calling Sequence:

First, the driver starts initialization. When the initialization of application installation and transmission handles has been successfully completed, the number of controllable cameras can be obtained.

When the application starts, it calls the camera initialization function to perform initialization. After the initialization function is successfully used, other functions can be called normally for execution.

The camera termination function is used for closing the program. It is the function to be executed when a camera is suspended, or the resources are released and no longer control the camera. For example, when exiting the application and the termination function is called, the calls of other performance functions will not be executed until after the initialization function is called again.

4.1.3 Driver Initialization:

Driver uses [TUCAM_Api_Init](#) function for the initialization operation. This function initialises

frame grabbers and controls the camera.

4.1.4 Camera Initialization:

Camera initialization uses `TUCAM_Dev_Open` function. This function retrieves the necessary camera handles to use the input parameters for other functions.

4.1.5 Camera Product Information:

When the camera is opened after calling `TUCAM_Dev_Open` function, the camera product information can be obtained through the camera handle.

```
// enumerate information id
typedef enum
{
    TUIDI_BUS           = 0x01,           // USB interface type: USB2.0/USB3.0
    TUIDI_VENDOR        = 0x02,           // Manufacturer ID
    TUIDI_PRODUCT        = 0x03,           // Product ID
    TUIDI_VERSION_API    = 0x04,           // TUCAM- API Version number
    TUIDI_VERSION_FRMW   = 0x05,           //Firmware version number
    TUIDI_VERSION_FPGA   = 0x06,           // FPGA Version number (reserved)
    TUIDI_VERSION_DRIVER = 0x07,           // Driver version number (reserved)
    TUIDI_TRANSFER_RATE  = 0x08,           // USB transfer rate
    TUIDI_CAMERA_MODEL   = 0x09,           // Camera model (string type)
    TUIDI_ENDINFO        = 0x0A,           // Product information ID end bit
}TUCAM_IDINFO;
```

Example

```
TUIDI_BUS           = 0x300,           // USB3.0
TUIDI_VENDOR        = 0x5453,          // TUCSEN
TUIDI_PRODUCT        = 0x6404,          // Dhyana 400D
TUIDI_VERSION_API    = "1.0.0.1",       // "1.0.0.1"
TUIDI_VERSION_FRMW   = -230244288 ,    // "f246c040"
TUIDI_VERSION_FPGA   = 0,               // the FPGA version (reserved)
TUIDI_VERSION_DRIVER = "1.2.3.10",     // "1.2.3.10"
TUIDI_TRANSFER_RATE  = 292,             // 292MB / Sec
TUIDI_CAMERA_MODEL   = "Dhyana 400D",  // "Dhyana 400D"
```

4.1.6 Termination Procedure:

The termination camera procedure uses `TUCAM_Dev_Close` function. Calling of this function releases the obtained port and the resources used for the camera frame. When this function is called, the camera will no longer be controlled.

4.1.7 Sample Code:

```

1.  int main (int argc, char** argv)
2.  {
3.      TUCAM_INIT  itApi; // Initializing SDK environment parameters
4.      TUCMA_OPEN  opCam; // Open camera parameters
5.
6.      itApi.pstrConfigPath = NULL;
7.      itApi.uiCamCount = 0;
8.      if (TUCAMRET_SUCCESS != TUCAM_Api_Init(&itApi))
9.      {
10.         // Initializing SDK API environment failed
11.         return 0;
12.      }
13.
14.      if (0 == itApi.uiCamCount)
15.      {
16.         // No camera
17.         return 0;
18.      }
19.
20.      opCam.hIdxTUCam = 0;
21.      opCam.uiIdxOpen = 0;
22.
23.      if (TUCAMRET_SUCCESS != TUCAM_Dev_Open(&opCam))
24.      {
25.         // Failed to open camera
26.         return 0;
27.      }
28.
29.      // Application can use the handle of opCam.hIdxTUCam
30.
31.      TUCAM_Dev_Close(opCam.hIdxTUCam); // Close camera
32.      TUCAM_Api_Uninit();                // Initializing SDK API environment
33.

```

```

34.     return 0;
35. }
36.

```

4.2 Performance Obtaining and Setting

4.2.1 Interface:

// Capability control see enumerate TUCAM_IDCAPA

TUCAMRET TUCAM_Capa_GetAttr (HDTUCAM hTUCam, PTUCAM_CAPA_ATTR pAttr);

TUCAMRET TUCAM_Capa_GetValue (HDTUCAM hTUCam, INT32 nCapa, INT32 *pnVal);

TUCAMRET TUCAM_Capa_SetValue (HDTUCAM hTUCam, INT32 nCapa, INT32 nVal);

TUCAMRET TUCAM_Capa_GetValueText (HDTUCAM hTUCam, PTUCAM_VALUE_TEXT pVal);

4.2.2 Calling Sequence:

The obtaining and setting are generally completed before or after the camera capture. If the function is set to be called when the data is captured, in some cases, it may return the error code TUCAMRET.

4.2.3 Performance Index:

// enumerate capability id

typedef enum

```

{
    TUIDC_RESOLUTION           = 0x00,           // Resolution
    TUIDC_PIXELCLOCK           = 0x01,           // Pixel clock
    TUIDC_BITOFDEPTH           = 0x02,           // Data bit wide
    TUIDC_ATEXPOSURE            = 0x03,           // Automatic exposure
    TUIDC_HORIZONTAL            = 0x04,           // Horizontal Flip
    TUIDC_VERTICAL              = 0x05,           // Vertical Flip
    TUIDC_ATWBBALANCE           = 0x06,           // Automatic white balance (color
camera)
    TUIDC_FAN_GEAR              = 0x07,           // Fan level (refrigeration camera)
    TUIDC_ENDCAPABILITY         = 0x08,           // Performance ID bit end
}TUCAM_IDCAPA;

```

Note: If the camera does not support the performance ID, it will return the error code TUCAMRET_NOT_SUPPORT.

4.2.3 Sample Code:

```

1. // Exampled with the resolution TUIDC_RESOLUTION
2. // Obtain the resolution scope
3. void GetResolutionRange()
4. {
5.     TUCAM_CAPA_ATTR attrCapa;
6.     TUCAM_VALUE_TEXT valText;
7.
8.     char szRes[64] = {0};
9.     valText.nTextSize = 64;
10.    valText.pText = &szRes[0];
11.    attrCapa.idCapa = TUIDC_RESOLUTION;
12.    if (TUCAMRET_SUCCESS == TUCAM_Capa_GetAttr(opCam.hIdxTUCam, &attrCapa))
13.    {
14.        // Obtain the number of resolution
15.        int nCnt = attrCapa.nValMax - attrCapa.nValMin + 1;
16.        valText.nID = TUIDC_RESOLUTION;
17.
18.        for (int i=0; i<nCnt; ++i)
19.        {
20.            valText.dbValue = i;
21.            TUCAM_Capa_GetValueText(opCam.hIdxTUCam, &valText);
22.            szRes = valText.pText;
23.            // Add resolution text to the drop-down menu
24.        }
25.    }
26. }
27.
28. // Obtain the current resolution
29. void GetCurrentResolution()
30. {
31.     int nVal = 0;
32.
33.     if (TUCAMRET_SUCCESS == TUCAM_Capa_GetValue(opCam.hIdxTUCam, \
34.                                                    TUIDC_RESOLUTION, \
35.                                                    &nVal))
36.     {
37.         // nVal Returns to the current resolution index
38.     }
39. }
40.
41. // Set the current resolution
42. void SetCurrentResolution(int nIdxRes)

```

```

43. {
44.     TUCAM_Capa_SetValue(opCam.hIdxTUCam, TUIDC_RESOLUTION, nIdxRes);
45. }
46.

```

4.3 Property Extraction and Settings

4.3.1 Interface:

// Property control see enumerate TUCAM_IDPROP

```

TUCAMRET TUCAM_Prop_GetAttr (HDTUCAM hTUCam, PTUCAM_PROP_ATTR pAttr);
TUCAMRET TUCAM_Prop_GetValue (HDTUCAM hTUCam, INT32 nProp, DOUBLE *pdbVal, INT32 nChn);
TUCAMRET TUCAM_Prop_SetValue (HDTUCAM hTUCam, INT32 nProp, DOUBLE dbVal, INT32
nChn);
TUCAMRET TUCAM_Prop_GetValueText (HDTUCAM hTUCam, PTUCAM_VALUE_TEXT pVal, INT32
nChn );

```

4.3.2 Calling Sequence:

The obtaining and setting are generally completed before or after the camera capture. If the function is set to be called when the data is captured, in some cases, it may return the error code TUCAMRET.

4.3.3 Property Index:

// enumerate property id

typedef enum

```

{
    TUIDP_GLOBALGAIN      = 0x00,          // Global gain
    TUIDP_EXPOSURETM      = 0x01,          // Exposure time
    TUIDP_BRIGHTNESS      = 0x02,          // Brightness (valid under AE status)
    TUIDP_BLACKLEVEL      = 0x03,          // Black level
    TUIDP_TEMPERATURE     = 0x04,          // Temperature
    TUIDP_SHARPNESS       = 0x05,          // Sharpness
    TUIDP_NOISELEVEL      = 0x06,          // Noise level
    TUIDP_HDR_KVALUE      = 0x07,          // HDR K value (for sCMOS cameras)

```

// image process property

```

    TUIDP_GAMMA           = 0x08,          // Gamma
    TUIDP_CONTRAST        = 0x09,          // Contrast
    TUIDP_LFTLEVELS      = 0x0A,          // Left Levels

```

```

TUIDP_RGTLEVELS      = 0x0B,          // Right Levels
TUIDP_CHNLGAIN        = 0x0C,          // Channel gain (for color cameras)
TUIDP_SATURATION      = 0x0D,          // Saturation (for color cameras)
TUIDP_ENDPROPERTY     = 0x0E,          // Property ID end bit
}TUCAM_IDPROP;

```

Note: If the camera does not support the Property ID, it will return the error code **TUCAMRET_NOT_SUPPORT**.

4.3.3 Sample Code:

```

1.  // Exemplified with exposure time
2.  // Obtain exposure time range
3.  void GetExposureTimeRange()
4.  {
5.      TUCAM_PROP_ATTR attrProp;
6.
7.      attrProp.nIdxChn = 0;    // Current channel
8.      attrProp.idProp = TUIDP_EXPOSURETM;
9.
10.     if (TUCAMRET_SUCCESS == TUCAM_Prop_GetAttr(opCam.hIdxTUCam, &attrProp))
11.     {
12.         // Exposure time range
13.         attrProp.dbValMin;    // Minimum exposure time
14.         attrProp.dbValMax;    // Maximum exposure time
15.         attrProp.dbValDft;    // Default exposure time
16.         attrProp.dbValStep;   // Exposure time step
17.     }
18. }
19.
20. // Obtain current exposure time
21. void GetCurrentExposureTime()
22. {
23.     double dbVal = 1.0f;
24.
25.     if (TUCAMRET_SUCCESS == TUCAM_Prop_GetValue(opCam.hIdxTUCam, \
26.                                                  TUIDP_EXPOSURETM, \
27.                                                  &dbVal))
28.     {
29.         // dbVal returns the current exposure time in ms
30.     }
31. }
32.

```

```

33. // Set current exposure time
34. void SetCurrentExposureTime(double dbVal)
35. {
36.     TUCAM_Prop_SetValue(opCam.hIdxTUCam, TUIDP_EXPOSURETM, dbVal);
37. }
38.

```

4.4 Memory management

4.4.1 Interface:

```

// Buffer control see structure TUCAM_FRAME
TUCAMRET TUCAM_Buf_Alloc (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
TUCAMRET TUCAM_Buf_Release (HDTUCAM hTUCam);
TUCAMRET TUCAM_Buf_AbortWait (HDTUCAM hTUCam);
TUCAMRET TUCAM_Buf_WaitForFrame (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
TUCAMRET TUCAM_Buf_CopyFrame (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);

```

4.4.2 Calling Sequence:

For memory allocation and release, the memory allocation **TUCAM_Buf_Alloc** must be called before **TUCAM_Cap_Start** starts data capture; while memory release **TUCAM_Buf_Release** must be called after **TUCAM_Cap_Stop** stops the data capture.

For the data capture, **TUCAM_Buf_WaitForFrame** must be called after **TUCAM_Cap_Start** starts the data capture

Wait for the completion of data capture, and data in different formats can be copied via **TUCAM_Buf_CopyFrame**.

If there are calls of data waiting and data copying, prior to stopping the data capture, **TUCAM_Buf_AbortWait** is called to stop data waiting, and then **TUCAM_Cap_Stop** is called to stop data capture.

4.4.3 Frame Structure:

```

// the camera frame structure
typedef struct _tagTUCAM_FRAME
{
    // TUCAM_Buf_WaitForFrame Using this structure enables a number of member variables of
    // different call directions.
    // "input" means that the application is called before setting

```

```

// "output" means the value required is returned after the program calls the interface

CHAR szSignature[8];    // [out] Copyright information

// The based information
USHORT usHeader;        // [out] Head size of frame
USHORT usOffset;        // [out] Head offset of frame data (generally in the same size as the
head)
USHORT usWidth;         // [out] Frame width
USHORT usHeight;        // [out] Frame height
UINT32 uiWidthStep;     // [out] Frame width step

UCHAR ucDepth;          // [out] Frame data bit depth
UCHAR ucFormat;         // [out] Frame data format
UCHAR ucChannels;       // [out] Number of frame data channel
UCHAR ucElemBytes;      // [out] Frame pixel bytes
UCHAR ucFormatGet;      // [in] Obtain frame format (refer to TUFrm_FORMATS)

UINT32 uiIndex;         // [in/out] Current frame number
UINT32 uiImgSize;       // [out] Frame size
UINT32 uiRsdSize;       // [in] Number of reserved frames (number of frames needed, trigger
use)
UINT32 uiHstSize;       // [out] Frame histogram frame size (reserved bit)

PUCHAR pBuffer;        // [in/out] Frame buffer
} TUCAM_FRAME, *PTUCAM_FRAME;

// enumerate frame format
typedef enum
{
    TUFrm_FMT_RAW          = 0x10,          // Data of Raw format
    TUFrm_FMT_USUAI        = 0x11,          // General data (8bit/16bit, monochrome,
color)
    TUFrm_FMT_RGB888       = 0x12,          // Data of RGB888 for display
}TUFrm_FORMATS;

```

4.4.4 Sample Code:

```

1. TUCAM_FRAME m_frame;    // Frame object
2. HANDLE m_hThdGrab;      // Event handler of picture grabbing
3. BOOL m_bLiving;         // Confirm to grab picture
4.

```

```

5.  BOOL CDlgTUCam::StartCapture()
6.  {
7.      m_frame.pBuffer      = NULL;
8.      m_frame.ucFormatGet = TUFRM_FMT_RGB888; // Frame data format (RGB888)
9.      m_frame.uiRsdSize    = 1; // Number of frames captured once
    (TUCCM_TRIGGER_STANDARD may be greater than 1)
10.
11.     if (TUCAMRET_SUCCESS != TUCAM_Buf_Alloc(m_opCam.hIdxTUCam, &m_frame))
12.     {
13.         return FALSE;
14.     }
15.
16.     if (TUCAMRET_SUCCESS != TUCAM_Cap_Start(m_opCam.hIdxCam,
    TUCCM_SEQUENCE))
17.     {
18.         TUCAM_Buf_Release(m_opCam.hIdxTUCam);
19.         return FALSE;
20.     }
21.
22.     m_bLiving = TRUE;
23.     m_hThdGrab = CreateEvent(NULL, TRUE, FALSE, NULL);
24.     _beginthread(GrabThread, 0, this);
25.
26.     return TRUE;
27. }
28.
29. Void __cdecl CDlgTUCam::GrabThread(LPVOID IPParam)
30. {
31.     CDlgTUCam *pTuCam = (CDlgTUCam *)IPParam;
32.
33.     While (pTUCam->m_bLiving)
34.     {
35.         pTUCam->m_frame.ucFormatGet = TUFRM_FMT_RGB888;
36.         if(TUCAMRET_SUCCESS ==
    TUCAM_Buf_WaitForFrame(pTUCam->m_opCam.hIdxTUCam,\
37.                                     &pTUCam->m_frame))
38.         {
39.             // pTUCam->m_frame.pBuffer Return the captured image data in the format of
    TUFRM_FMT_RGB88
40.             // The data can be used to display
41.
42.             // Data in other formats can be obtained
43.             pTUCam->m_frame.ucFormatGet = TUFRM_FMT_USUAL;
44.

```

```

    if(TUCAMRET_SUCCESS==TUCAM_Buf_CopyFrame(pTUCam->m_opCam.hIdxTUCam,\
45.                                     &pTUCam->m_frame))
46.     {
47.         // pTUCam->m_frame.pBuffer Return the captured image data
48.     }
49. }
50. }
51.
52. SetEvent(pTUCam->m_hThdGrab);
53. _endthread();
54. }
55.
56. void CDlgTUCam::StopCapture()
57. {
58.     m_bLiving = FALSE;
59.     TUCAM_BUF_AbortWait();    // If calls TUCAM_Buf_WaitForFrame Interface
60.
61.     WaitForSingleObject(m_hThdGrab, INFINITE);    // Wait for exiting of picture grabbing
62.     CloseHandle(m_hThdGrab);
63.     m_hThdGrab = NULL;
64.
65.     TUCAM_Cap_Stop(m_opCam.hIdxTUCam);    // Stop data capture
66.     TUCAM_Buf_Release(m_opCam.hIdxTUCam);    // Release the allocated memory
67. }
68.

```

4.5 Capture Control

4.5.1 Interface:

// Capturing control

// ROI

TUCAMRET TUCAM_Cap_SetROI (HDTUCAM hTUCam, TUCAM_ROI_ATTR roiAttr);

TUCAMRET TUCAM_Cap_GetROI (HDTUCAM hTUCam, PTUCAM_ROI_ATTR pRoiAttr);

// Trigger

TUCAMRET TUCAM_Cap_SetTrigger (HDTUCAM hTUCam, TUCAM_TRIGGER_ATTR tgrAttr);

TUCAMRET TUCAM_Cap_GetTrigger (HDTUCAM hTUCam, PTUCAM_TRIGGER_ATTR pTgrAttr);

TUCAMRET TUCAM_Cap_DoSoftwareTrigger(HDTUCAM hTUCam); // in trigger mode

// Capturing

// uiMode see enumerate TUCAM_CAPTURE_MODES

TUCAMRET TUCAM_Cap_Start(HDTUCAM hTUCam, UINT32 uiMode);

TUCAMRET TUCAM_Cap_Stop (HDTUCAM hTUCam);

4.5.2 Calling Sequence:

Set ROI property [TUCAM_Cap_SetROI](#) and trigger property [TUCAM_Cap_SetTrigger](#), which needs to be called before start capturing data; if called when capturing data, the error code TUCAMRET may be returned.

4.5.2 Capture Mode Index:

```
// enumerate the capture mode
typedef enum
{
    TUCCM_SEQUENCE           = 0x00,           // Sequence mode for data capture
    TUCCM_TRIGGER_STANDARD   = 0x01,           // Standard trigger mode for data capture
    TUCCM_TRIGGER_SYNCHRONOUS = 0x02,           // Synchronization trigger mode for data
capture
    TUCCM_TRIGGER_GLOBAL     = 0x03,           // Global trigger mode for data capture
    TUCCM_TRIGGER_SOFTWARE   = 0x04,           // Software trigger mode for data
capture
}TUCAM_CAPTURE_MODES;
```

4.5.3 Sample Code:

```
1. // Set ROI mode
2. void SetROIMode()
3. {
4.     TUCAM_ROI_ATTR roiAttr;
5.     roiAttr.bEnable = TRUE;
6.     roiAttr.nVOffset= 100;
7.     roiAttr.nHOffset = 100;
8.     roiAttr.nWidth   = 800;
9.     roiAttr.nHeight  = 600;
10.
11.     TUCAM_Cap_SetROI(m_opCam.hIdxTUCam, roiAttr);
12.     TUCAM_Cap_Start(m_opCam.hIdxCam, TUCCM_SEQUENCE); // Sequence mode
        (stream mode)
13.
14.     // Refer to memory management sample code for data obtaining
15. }
16.
17. // Set trigger mode
```

```

18. void SetTriggerMode()
19. {
20.     TUCAM_TRIGGER_ATTR tgrAttr;
21.
22.     tgrAttr.nTgrMode = TUCCM_TRIGGER_STANDARD;    // Standard trigger mode
23.     tgrAttr.nExpMode = TUCTE_EXPTM;               // Exposure mode
24.     tgrAttr.nEdgeMode= TUCTE_RISING;              // Stimulate rising edge
25.     tgrAttr.nFrames  = 1;                         // Trigger one frame
26.     tgrAttr.nDelayTm = 0;                         // Delay 0 ms
27.
28.     TUCAM_Cap_SetTrigger(m_opCam.hIdxTUCam, tgrAttr);
29.     TUCAM_Cap_Start(m_opCam.hIdxCam, TUCCM_STANDARD);    // Standard trigger
mode
30.
31.     // Refer to memory management sample code for data obtaining
32. }
33.

```

4.6 File Control

4.6.1 Interface:

```

// File control
// Image
TUCAMRET  TUCAM_File_SaveImage (HDTUCAM hTUCam, TUCAM_FILE_SAVE fileSave);
// Video
TUCAMRET  TUCAM_Rec_Start(HDTUCAM hTUCam, TUCAM_REC_SAVE recSave);
TUCAMRET  TUCAM_Rec_AppendFrame(HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
TUCAMRET  TUCAM_Rec_Stop (HDTUCAM hTUCam);

```

4.6.2 Calling Sequence:

Starting of video recording **TUCAM_Rec_Start** needs to be called after **TUCAM_Cap_Start** starts capturing data, and the set capture mode must be **TUCCM_SEQUENCE** mode. In the video recording process, by calling **TUCAM_Rec_AppendFrame**, the image data are written into the file, the end of the video recording calls **TUCAM_Rec_Stop** to end the video recording process.

4.6.3 File Structure:

```

// the file save structure

```

```

typedef struct _tagTUCAM_FILE_SAVE
{
    INT32    nSaveFmt;           // [in] Format of saved file (refer to TUIMG_FORMATS)
    PCHAR    pstrSavePath;       // [in] Path of saved file (including file name, but not including
the extension name)
    PTUCAM_FRAME pFrame;         // [in] Frame structure pointer
} TUCAM_FILE_SAVE, *PTUCAM_FILE_SAVE;

// the record save structure
typedef struct _tagTUCAM_REC_SAVE
{
    INT32    nCodec;             // [in] Codec type
    PCHAR    pstrSavePath;       // [in] Path of saved file (including file name, but not including
the extension name)
    Float    fFps;               // [in] Current frame rate (video frame rate)
} TUCAM_REC_SAVE, *PTUCAM_REC_SAVE;

```

4.6.4 Sample Code:

```

1.  // Save image file
2.  void SaveImage()
3.  {
4.      m_frame.ucFormatGet = TUFRTM_FMT_USUAL;
5.      if(TUCAMRET_SUCCESS==TUCAM_Buf_WaitForFrame(m_opCam.hIdxTUCam,
&m_frame))
6.      {
7.          TUCAM_FILE_SAVE fileSave;
8.          fileSave.nSaveFmt    = TUFMT_TIF;        // Save Tiff format
9.          fileSave.pFrame      = &m_frame;         // Frame pointer needs to be saved
10.         fileSave.pstrSavePath = "C:\\image";      // Path including file name (not including
extension name)
11.
12.         if (TUCAMRET_SUCCESS == TUCAM_File_SaveImage(m_opCam.hIdxTUCam,
fileSave))
13.         {
14.             // Image file saved successfully
15.         }
16.     }
17. }
18.
19. // Save video file
20. void StartRecording()
21. {

```



```

22.   TUCAM_REC_SAVE   recSave;
23.   recSave.fFps      = 15.0f;           // Frame rate needs to be saved
24.   recSave.nCodec    = m_dwFccHandler;
25.   recSave.pstrSavePath = "C:\\TUVVideo.avi" // Full path
26.
27.   if (TUCAMRET_SUCCESS == TUCAM_Rec_Start(m_opCam.hIdxTUCam, recSave))
28.   {
29.       // Start video recording. . .
30.   }
31. }
32.
33. Void AppendFrame()
34. {
35.     m_frame.ucFormatGet = TUFRM_FMT_RGB888;
36.     if(TUCAMRET_SUCCESS==TUCAM_Buf_WaitForFrame(m_opCam.hIdxTUCam,
37.         &m_frame))
38.     {
39.         TUCAM_Rec_AppendFrame(m_opCam.hIdxTUCam, &m_frame);
40.     }
41. }
42. void StopRecording()
43. {
44.     TUCAM_Rec_Stop(m_opCam.hIdxTUCam);
45. }
46.

```

4.7 Extending Control

4.7.1 Interface:

// Extended control

TUCAMRET TUCAM_Reg_Read (HDTUCAM hTUCam, TUCAM_REG_RW regRW);

TUCAMRET TUCAM_Reg_Write(HDTUCAM hTUCam, TUCAM_REG_RW regRW);

4.7.2 Calling Sequence:

Read/write register of **TUCAM_Reg_Read** / **TUCAM_Reg_Write** must be called after **TUCAM_Dev_Open** opens the camera; if the register cannot be read and written after **TUCAM_Dev_Close** turns off the camera, the camera must be reopened.

4.7.3 File Structure:

```
// the register read/write struct
typedef struct _tagTUCAM_REG_RW
{
    INT32    nRegFmt;           // [in] Type of read/write register (refer to TUREG_TYPE)
    PCHAR    pBuf;             // [in/out] Pointer to the buffer
    INT32    nBufSize;         // [in] Buffer size
} TUCAM_REG_RW, *PTUCAM_REG_RW;
```

4.7.4 Sample Code:

```
1.  // Read register data
2.  void ReadRegisterData()
3.  {
4.      char cSN[TUSN_SIZE] = {0};
5.      TUCAM_REG_RW regRW;
6.
7.      regRW.nRegType = TUREG_SN;
8.      regRW.pBuf      = &cSN[0];
9.      regRW.nBufSize  = TUSN_SIZE;
10.
11.     if (TUCAMRET_SUCCESS == TUCAM_Reg_Read(m_opCam.hIdxTUcam, regRW))
12.     {
13.         // Obtain SN data
14.     }
15. }
16.
17. // Write register data
18. void WriteRegisterData()
19. {
20.     char cSN[TUSN_SIZE] = {'S', 'N', '1', '2', '3', '4', '5', '6'}; // "SN123456"
21.     TUCAM_REG_RW regRW;
22.
23.     regRW.nRegType = TUREG_SN;
24.     regRW.pBuf      = &cSN[0];
25.     regRW.nBufSize  = TUSN_SIZE;
26.
27.     if (TUCAMRET_SUCCESS == TUCAM_Reg_Write(m_opCam.hIdxTUcam, regRW))
28.     {
29.         // SN successfully written into register
30.     }
```

```
31. }
```

5. Reference

5.1 Types and Constants

5.1.1 TUCAMRET Error code:

TUCAMRET_SUCCESS	= 0x00000001,	// No error, general success code
TUCAMRET_FAILURE	= 0x80000000,	// Error
// initialization error		
TUCAMRET_NO_MEMORY	= 0x80000101,	// Not enough memory
TUCAMRET_NO_RESOURCE	= 0x80000102,	// Not enough resources (not including memory)
TUCAMRET_NO_MODULE	= 0x80000103,	// No sub-module
TUCAMRET_NO_DRIVER	= 0x80000104,	// No driver
TUCAMRET_NO_CAMERA	= 0x80000105,	// No camera
TUCAMRET_NO_GRABBER	= 0x80000106,	// No picture grabber
TUCAMRET_NO_PROPERTY	= 0x80000107,	// No alternative property ID
TUCAMRET_FAILOPEN_CAMERA	= 0x80000110,	// Failed to open camera
TUCAMRET_FAILOPEN_BULKIN	= 0x80000111,	// Failed to open batch transmission input terminal
TUCAMRET_FAILOPEN_BULKOUT	= 0x80000112,	// Failed to open batch transmission output terminal
TUCAMRET_FAILOPEN_CONTROL	= 0x80000113,	// Failed to open control endpoint
TUCAMRET_FAILCLOSE_CAMERA	= 0x80000114,	// Failed to close camera
TUCAMRET_FAILOPEN_FILE	= 0x80000115,	// Failed to open file
// status error		
TUCAMRET_INIT	= 0x80000201,	// API needs to initialize state.
TUCAMRET_BUSY	= 0x80000202,	// API busy
TUCAMRET_NOT_INIT	= 0x80000203,	// API is not initialized
TUCAMRET_EXCLUDED	= 0x80000204,	// Some resources are used exclusively
TUCAMRET_NOT_BUSY	= 0x80000205,	// API is not busy
TUCAMRET_NOT_READY	= 0x80000206,	// API is not standby
// wait error		
TUCAMRET_ABORT	= 0x80000207,	// Processing aborted
TUCAMRET_TIMEOUT	= 0x80000208,	// Timeout
TUCAMRET_LOSTFRAME	= 0x80000209,	// Frame loss
TUCAMRET_MISSFRAME	= 0x8000020A,	// Frame loss but it is caused by underlying driver issue
// calling error		

TUCAMRET_INVALID_CAMERA	= 0x80000301,	// Invalid camera
TUCAMRET_INVALID_HANDLE	= 0x80000302,	// Invalid camera handle
TUCAMRET_INVALID_OPTION	= 0x80000303,	// Invalid configuration value
TUCAMRET_INVALID_IDPROP	= 0x80000304,	// Invalid property ID
TUCAMRET_INVALID_IDCAPA	= 0x80000305,	// Invalid performance ID
TUCAMRET_INVALID_IDPARAM	= 0x80000306,	// Invalid parameter ID
TUCAMRET_INVALID_PARAM	= 0x80000307,	// Invalid parameter
TUCAMRET_INVALID_FRAMEIDX	= 0x80000308,	// Invalid frame IDX
TUCAMRET_INVALID_VALUE	= 0x80000309,	// Invalid value
TUCAMRET_INVALID_EQUAL	= 0x8000030A,	// Values are equal, but the parameter is invalid
TUCAMRET_INVALID_CHANNEL	= 0x8000030B,	// Specified channel of property ID, but the channel is not valid
TUCAMRET_INVALID_SUBARRAY	= 0x8000030C,	// Value of the sub-array is invalid
TUCAMRET_INVALID_VIEW	= 0x8000030D,	// Invalid display window handle
TUCAMRET_INVALID_PATH	= 0x8000030E,	// Invalid file path
TUCAMRET_NO_VALUETEXT	= 0x80000310,	// Text without property value
TUCAMRET_OUT_OF_RANGE	= 0x80000311,	// Value is out of range
TUCAMRET_NOT_SUPPORT	= 0x80000312,	// Unsupported features or properties
TUCAMRET_NOT_WRITABLE	= 0x80000313,	// Property unwritable
TUCAMRET_NOT_READABLE	= 0x80000314,	// Property unreadable
TUCAMRET_WRONG_HANDSHAKE	= 0x80000410,	// Error occurred while retrieving the error code
TUCAMRET_NEWAPI_REQUIRED	= 0x80000411,	// Old API does not support, only new API supports
TUCAMRET_ACCESSDENY	= 0x80000412,	// Camera cannot be accessed under certain state
TUCAMRET_NO_CORRECTIONDATA	= 0x80000501,	// No color dot correction of data.
// camera or bus trouble		
TUCAMRET_FAIL_READ_CAMERA	= 0x83001001,	// Failed to read from camera
TUCAMRET_FAIL_WRITE_CAMERA	= 0x83001002,	// Failed to write into camera
TUCAMRET_OPTICS_UNPLUGGED	= 0x83001003,	// Unplugged

5.1.2 TUCAM_IDINFO Product Information Code:

TUIDI_BUS	= 0x01,	// USB interface type: USB2.0/USB3.0
TUIDI_VENDOR	= 0x02,	// Manufacturer ID

TUIDI_PRODUCT	= 0x03,	// Product ID
TUIDI_VERSION_API	= 0x04,	// TUCAM- API Version number
TUIDI_VERSION_FRMW	= 0x05,	// Firmware version number
TUIDI_VERSION_FPGA	= 0x06,	// FPGA Version number (reserved)
TUIDI_VERSION_DRIVER	= 0x07,	// Driver version number (reserved)
TUIDI_TRANSFER_RATE	= 0x08,	// USB transfer rate
TUIDI_CAMERA_MODEL	= 0x09,	// Camera model (string type)
TUIDI_ENDINFO	= 0x0A,	// Product information ID end bit

5.1.3 TUCAM_IDCAPA Performance Code:

TUIDC_RESOLUTION	= 0x00,	// Resolution
TUIDC_PIXELCLOCK	= 0x01,	// Pixel clock
TUIDC_BITOFDEPTH	= 0x02,	// Data bit wide
TUIDC_ATEXPOSURE	= 0x03,	// Automatic exposure
TUIDC_HORIZONTAL	= 0x04,	// Horizontal Flip
TUIDC_VERTICAL	= 0x05,	// Vertical Flip
TUIDC_ATWBBALANCE	= 0x06,	// Automatic white balance (color camera)
TUIDC_FAN_GEAR	= 0x07,	// Fan level (refrigeration camera)
TUIDC_ENDCAPABILITY	= 0x08,	// Performance ID bit end

5.1.4 TUCAM_IDPROP Property Code:

TUIDP_GLOBALGAIN	= 0x00,	// Global gain
TUIDP_EXPOSURETM	= 0x01,	// Exposure time
TUIDP_BRIGHTNESS	= 0x02,	// Brightness (valid under AE status)
TUIDP_BLACKLEVEL	= 0x03,	// Black level
TUIDP_TEMPERATURE	= 0x04,	// Temperature
TUIDP_SHARPNESS	= 0x05,	// Sharpness
TUIDP_NOISELEVEL	= 0x06,	// Noise level
TUIDP_HDR_KVALUE	= 0x07,	// HDR K value (for sCMOS cameras)
// image process property		
TUIDP_GAMMA	= 0x08,	// Gamma
TUIDP_CONTRAST	= 0x09,	// Contrast
TUIDP_LFTLEVELS	= 0x0A,	// Left Levels
TUIDP_RGTLEVELS	= 0x0B,	// Right Levels
TUIDP_CHNLGAIN	= 0x0C,	// Channel gain (for color cameras)
TUIDP_SATURATION	= 0x0D,	// Saturation (for color cameras)
TUIDP_ENDPROPERTY	= 0x0E,	// Property ID end bit

5.1.5 TUCAM_CAPTURE_MODES Capture Mode Code:

TUCCM_SEQUENCE	= 0x00,	// Sequence mode (stream mode)
TUCCM_TRIGGER_STANDARD	= 0x01,	// Standard trigger mode
TUCCM_TRIGGER_SYNCHRONOUS	= 0x02,	// Synchronization trigger mode
TUCCM_TRIGGER_GLOBAL	= 0x03,	// Global trigger mode
TUCCM_TRIGGER_SOFTWARE	= 0x04,	// Software trigger mode

5.1.6 TUIMG_FORMATS Image Format Code:

TUFMT_RAW	= 0x01,	// RAW format
TUFMT_TIF	= 0x02,	// TIFF format
TUFMT_PNG	= 0x04,	// PNG format
TUFMT_JPG	= 0x08,	// JPEG format
TUFMT_BMP	= 0x10,	// BMP format

5.1.7 TUREG_TYPE Register Type Code:

TUREG_SN	= 0x01,	// Read and write SN code register
TUREG_DATA	= 0x02,	// Read and write DATA register
(reserved)		

5.1.8 TUCAM_TRIGGER_EXP Trigger Exposure Mode Code:

TUCTE_EXPTM	= 0x00,	// Exposure time mode for trigger
TUCTE_WIDTH	= 0x01,	// Level width mode for trigger

5.1.9 TUCAM_TRIGGER_EDGE Trigger Excitation Edge Code:

TUCTD_RISING	= 0x01,	// Stimulate rising edge
TUCTD_FALLING	= 0x00,	// Stimulate falling edge

5.1.10 TUFRM_FORMATS Frame Format Code:

TUFRM_FMT_RAW	= 0x10,	// RAW data format
TUFRM_FMT_USUI	= 0x11,	// Common format data (8bit/16bit, monochrome/color)
TUFRM_FMT_RGB888	= 0x12,	// RGB888 format data (can be used for display)

5.2 Structure

5.2.1 Initialization:

```
// the camera initialize structure
typedef struct _tagTUCAM_INIT
{
    UINT32    uiCamCount;           // [out] Return the number of the current
    connected cameras
    PCHAR     pstrConfigPath;       // [in] Camera parameters input saving path
}TUCAM_INIT, *PTUCAM_INIT;
```

5.2.2 Open Camera:

```
// the camera open structure
typedef struct _tagTUCAM_OPEN
{
    UINT32    uIdxOpen;             // [in] Input serial number of camera to be
    opened
    HDTUCAM   hIdxTUCam;           // [out] Output handles of opened cameras
}TUCAM_OPEN, *PTUCAM_OPEN;
```

5.2.3 Camera Information:

```
// the camera value text structure
typedef struct _tagTUCAM_VALUE_INFO
{
    INT32     nID;                  // [in] Information ID TUCAM_IDINFO
    INT32     nValue;              // [in] Information value
    PCHAR     pText;               // [in/out] Pointer pointing to text data
    INT32     nTextSize;           // [in] Text buffer size
}TUCAM_VALUE_INFO, *PTUCAM_VALUE_INFO;
```

5.2.4 Performance/Property Value Text:

```
// the camera value text structure
typedef struct _tagTUCAM_VALUE_TEXT
{
    INT32     nID;                  // [in] ID TUCAM_IDPROP / TUCAM_IDCAPA
    DOUBLE    dbValue;             // [in] Performance/property value
    PCHAR     pText;               // [in/out] Pointer pointing to text data
    INT32     nTextSize;           // [in] Text buffer size
}TUCAM_VALUE_TEXT, *PTUCAM_VALUE_TEXT;
```


5.2.5 Performance Property:

// the camera capability attribute

```
typedef struct _tagTUCAM_CAPA_ATTR
{
    INT32    idCapa;                // [in] ID TUCAM_IDCAPA

    INT32    nValMin;              // [out] Minimum value
    INT32    nValMax;              // [out] Maximum value
    INT32    nValDft;              // [out] Default value
    INT32    nValStep;             // [out] Step length
}TUCAM_CAPA_ATTR, *PTUCAM_CAPA_ATTR;
```

5.2.6 Property Attribute:

// the camera property attribute

```
typedef struct _tagTUCAM_PROP_ATTR
{
    INT32    idProp;               // [in] ID TUCAM_IDPROP
    INT32    nIdxChn;              // [in/out] Index number of current channel
    DOUBLE   dbValMin;             // [out] Minimum value
    DOUBLE   dbValMax;             // [out] Maximum value
    DOUBLE   dbValDft;             // [out] Default value
    DOUBLE   dbValStep;            // [out] Step length
}TUCAM_PROP_ATTR, *PTUCAM_PROP_ATTR;
```

5.2.7 ROI Attribute:

// the camera ROI attribute

```
typedef struct _tagTUCAM_ROI_ATTR
{
    BOOL     bEnable;              // [in/out] ROI enable

    INT32    nHOffset;             // [in/out] Horizontal offset
    INT32    nVOffset;             // [in/out] Vertical offset
    INT32    nWidth;               // [in/out] ROI width
    INT32    nHeight;              // [in/out] ROI height
}TUCAM_ROI_ATTR, *PTUCAM_ROI_ATTR;
```

5.2.8 Trigger Attribute:

// the camera trigger attribute

```
typedef struct _tagTUCAM_TRIGGER_ATTR
{

```

```

INT32  nTgrMode;           // [in/out] Trigger mode
INT32  nExpMode;           // [in/out] Exposure mode value [0,1] 0: Exposure
Time   1: level width
INT32  nEdgeMode;          // [in/out] Edge excitation mode [0, 1] 0: Falling
edge   1: Rising edge
INT32  nDelayTm;           // [in/out] Trigger delay time ms
INT32  nFrames;            // [in/out] Number of output frames of one trigger
}TUCAM_TRIGGER_ATTR, *PTUCAM_TRIGGER_ATTR;

```

5.2.9 Frame Structure:

```

// the camera frame structure
typedef struct _tagTUCAM_FRAME
{
    CHAR szSignature[8];    // [out] Copyright information

    // The based information
    USHORT usHeader;        // [out] Head size of frame
    USHORT usOffset;        // [out] Offset size of frame data
    USHORT usWidth;         // [out] Width of frame image
    USHORT usHeight;        // [out] Height of frame image
    UINT32 uiWidthStep;     // [out] Frame image width step

    UCHAR ucDepth;          // [out] Frame image data bit depth
    UCHAR ucFormat;         // [out] Frame image data format
    UCHAR ucChannels;       // [out] Number of frame image channel
    UCHAR ucElemBytes;      // [out] Frame image data byte
    UCHAR ucFormatGet;      // [in/out] Image format needs to be obtained
}TUCAM_FRAME, *PTUCAM_FRAME;

// TUCAM_FRAME_FORMATS

    UINT32 uiIndex;         // [out] Frame image serial number (reserved)
    UINT32 uiImgSize;       // [out] Frame image data size
    UINT32 uiRsdSize;       // [in] Number of frame needs to be obtained
    UINT32 uiHstSize;       // [out] Reserved field of frame image

    PUCAM_FRAME pBuffer;    // [in/out] Buffer pointing to frame data
} TUCAM_FRAME, *PTUCAM_FRAME;

```

5.2.10 File Saving:

```

// the file save structure
typedef struct _tagTUCAM_FILE_SAVE
{
    INT32  nSaveFmt;        // [in] Format of saved file refer to

```

TUIMG_FORMATS

```
PCHAR    pstrSavePath;           // [in] Path of saving (Not including the extension
name)
```

```
PTUCAM_FRAME pFrame;           // [in] Structure pointing to the frame
} TUCAM_FILE_SAVE, *PTUCAM_FILE_SAVE;
```

5.2.11 Video Recording Saving:

```
// the record save structure
```

```
typedef struct _tagTUCAM_REC_SAVE
```

```
{
```

```
    INT32    nCodec;           // [in] Codec type
```

```
    PCHAR    pstrSavePath;     // [in] Path of saved file including file name
```

```
    float    fFps;            // [in] Frame rate needs to be saved
```

```
} TUCAM_REC_SAVE, *PTUCAM_REC_SAVE;
```

5.2.12 Read/Write of Register

```
// the register read/write structure
```

```
typedef struct _tagTUCAM_REG_RW
```

```
{
```

```
    INT32    nRegType;         // [in] Type of read/write register refer to
```

```
TUREG_TYPE
```

```
    PCHAR    pBuf;             // [in/out] Point to the buffer of read/write contents
```

```
    INT32    nBufSize;         // [in] Buffer size
```

```
} TUCAM_REG_RW, *PTUCAM_REG_RW;
```

5.3 Functions

TUCAM_Api_Init

Description

Initialization of TUCAM-API library includes binding of driver and initialization of some internal resources, which is used before calling the other interfaces. The whole program only needs to call once.

Statement

```
TUCAMRET  TUCAM_Api_Init(PTUCAM_INIT pInitParam);
```

Parameters

PTUCAM_INIT pInitParam	Initialization of structure pointer, refer to structure TUCAM_INIT
------------------------	--

Error codes

TUCAMRET_INIT	TUCAM-API has been initialized
---------------	--------------------------------

Related Interfaces

TUCAM_Api_Uninit

TUCAM_Api_Uninit

Description

Uninstallation of TUCAM-API library includes the release of driver binding and some internal resources. It will be called once when the entire program is ended.

Statement

```
TUCAMRET  TUCAM_Api_Uninit ();
```

Parameters

No parameters input

Error codes

TUCAMRET_NOT_INIT

TUCAM-API not initialized

Related Interfaces

TUCAM_Api_Uninit

TUCAM_Dev_Open**Description**

Open the camera, the camera is in work mode after the call, which can respond to the calls of other interfaces. The camera should be ensured prior to that, that is, it should be after the initialization of calling [TUCAM_Api_Init](#).

Statement

```
TUCAMRET  TUCAM_Dev_Open (PTUCAM_OPEN pOpenParam);
```

Parameters

PTUCAM_OPEN pOpenParam Open camera structure pointer, refer to structure TUCAM_OPEN

Error Codes

TUCAMRET_NOT_INIT

TUCAM-API not initialized

TUCAMRET_INVALID_PARAM

Invalid parameter, when pOpenParam pointer is empty

TUCAMRET_OUT_OF_RANGE

Out of range, when the camera index needs to be opened exceeds the range of connected cameras

TUCAMRET_FAILOPEN_CAMERA

Failed to open camera

TUCAMRET_INVALID_CAMERA

Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit

TUCAM_Dev_Close

TUCAM_Dev_Close**Description**

Close the camera and the camera is in standby mode after called, and does not respond to calls from other interfaces.

Statement

```
TUCAMRET  TUCAM_Dev_Close ();
```

Parameters

No parameters input

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
-------------------	---------------------------

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open

TUCAM_Dev_GetInfo**Description**

Obtain the related information of camera, such as a USB port type, camera product number, API version number, firmware version number, camera type, etc. The camera should be ensured prior to that, and make sure the camera is opened, that is, it should be done after calling [TUCAM_Api_Init](#) initialization and [TUCAM_Api_Open](#).

Statement

```
TUCAMRET  TUCAM_Dev_GetInfo (HDTUCAM hTUCam, PTUCAM_VALUE_INFO pInfo);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_VALUE_INFO pInfo	Structure pointer of camera information value, refer to TUCAM_VALUE_INFO

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_PARAM	Invalid parameter, when product information code does not exist, refer to TUCAM_IDINFO
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open, TUCAM_Dev_Close

TUCAM_Capa_GetAttr

Description

Obtain the attribute value of performance parameters. The attributes obtained include the minimum value, maximum value, the default value and step of the parameter. The specific supported performance type can be referred to TUCAM_IDCAPA.

Statement

```
TUCAMRET  TUCAM_Capa_GetAttr (HDTUCAM hTUCam, PTUCAM_CAPA_ATTR pAttr);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_CAPA_ATTR pAttr	Camera performance attribute structure pointer, refer to TUCAM_CAPA_ATTR

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_IDCAPA	Invalid parameter, when product information code does not exist, refer to TUCAM_IDCAPA
TUCAMRET_INVALID_VALUE	Invalid value, when pAttr pointer is empty
TUCAMRET_NOT_SUPPORT	When the underlying request is not supported

TUCAMRET_INVALID_CAMERA Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit

TUCAM_Dev_Open, TUCAM_Dev_Close

TUCAM_Capa_GetValue, TUCAM_Capa_SetValue, TUCAM_Capa_GetValueText

TUCAM_Capa_GetValue

Description

Obtain the current attribute value of performance parameters. The specific supported performance type can be referred to TUCAM_IDCAPA.

Statement

```
TUCAMRET  TUCAM_Capa_GetValue (HDTUCAM hTUCam, INT32 nCapa, INT32 *pnVal);
```

Parameters

HDTUCAM hTUCam	Camera handle
INT32 nCapa	Camera performance attribute ID, refer to TUCAM_IDCAPA
INT32 *pnVal	Return the current value

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_IDCAPA	Invalid parameter, when product information code does not exist, refer to TUCAM_IDCAPA
TUCAMRET_INVALID_VALUE	Invalid value, when pAttr pointer is empty
TUCAMRET_NOT_SUPPORT	When the underlying request is not supported
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit

TUCAM_Dev_Open, TUCAM_Dev_Close

TUCAM_Capa_GetAttr, TUCAM_Capa_SetValue, TUCAM_Capa_GetValueText

TUCAM_Capa_SetValue

Description

Set the current attribute value of performance parameters. The specific supported performance type can be referred to TUCAM_IDCAPA.

Statement

```
TUCAMRET  TUCAM_Capa_SetValue (HDTUCAM hTUCam, INT32 nCapa, INT32 nVal);
```

Parameters

HDTUCAM hTUCam	Camera handle
INT32 nCapa	Camera performance attribute ID, refer to TUCAM_IDCAPA
INT32 nVal	Current value needs to be set

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_IDCAPA	Invalid parameter, when product information code does not exist, refer to TUCAM_IDCAPA
TUCAMRET_INVALID_VALUE	Invalid value, when pAttr pointer is empty
TUCAMRET_NOT_SUPPORT	When the underlying request is not supported
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Capa_GetAttr, TUCAM_Capa_GetValue, TUCAM_Capa_GetValueText

TUCAM_Capa_GetValueText

Description

Obtain the text information of current attribute value of performance parameters. The specific supported performance type can be referred to TUCAM_IDCAPA.

Statement

```
TUCAMRET  TUCAM_Capa_GetValueText (HDTUCAM hTUCam, PTUCAM_VALUE_TEXT pVal);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_VALUE_TEXT pVal	Obtain the text information structure pointer of performance parameters, TUCAM_VALUE_TEXT

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_FAILURE	Text buffer size is 0 or when pText pointer is empty
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Capa_GetAttr, TUCAM_Capa_SetValue, TUCAM_Capa_GetValue

TUCAM_Prop_GetAttr**Description**

Obtain the attribute value of attribute parameters. The attributes obtained include the minimum value, maximum value, the default value and step of the parameter. The specific supported performance type can be referred to TUCAM_IDPROP.

Statement

```
TUCAMRET  TUCAM_Prop_GetAttr (HDTUCAM hTUCam, PTUCAM_PROP_ATTR pAttr);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_PROP_ATTR pAttr	Camera performance attribute structure pointer, refer to TUCAM_PROP_ATTR

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
-------------------	---------------------------

TUCAMRET_INVALID_IDPROP	Invalid parameter, when product information code does not exist, refer to TUCAM_IDPROP
TUCAMRET_INVALID_VALUE	Invalid value, when pAttr pointer is empty
TUCAMRET_NOT_SUPPORT	When the underlying request is not supported
TUCAMRET_OUT_OF_RANGE	When the channel needs to be obtained exceeds the range
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit

TUCAM_Dev_Open, TUCAM_Dev_Close

TUCAM_Prop_GetValue, TUCAM_Prop_SetValue, TUCAM_Prop_GetValueText

TUCAM_Prop_GetValue

Description

Obtain the current attribute value of attribute parameters. The specific supported performance type can be referred to TUCAM_IDPROP.

Statement

```
TUCAMRET  TUCAM_Prop_GetValue (HDTUCAM hTUCam, INT32 nProp, DOUBLE *pdbVal, INT32
nChn = 0);
```

Parameters

HDTUCAM hTUCam	Camera handle
INT32 nProp	Camera performance attribute ID, refer to TUCAM_IDPROP
DOUBLE *pdbVal	Return the current value
INT32 nChn	Current channel needs to be obtained (the default is 0, monochrome camera is 0)

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_IDPROP	Invalid parameter, when product information code does not exist, refer to TUCAM_IDPROP
TUCAMRET_INVALID_VALUE	Invalid value, when pAttr pointer is empty
TUCAMRET_NOT_SUPPORT	When the underlying request is not supported
TUCAMRET_OUT_OF_RANGE	When the channel needs to be obtained exceeds the range
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Prop_GetAttr, TUCAM_Prop_SetValue, TUCAM_Prop_GetValueText

TUCAM_Prop_SetValue

Description

Set the current attribute value of attribute parameters. The specific supported performance type can be referred to TUCAM_IDPROP.

Statement

```
TUCAMRET  TUCAM_Prop_SetValue (HDTUCAM hTUCam, INT32 nProp, DOUBLE dbVal, INT32
nChn = 0);
```

Parameters

HDTUCAM hTUCam	Camera handle
INT32 nProp	Camera attribute ID, refer to TUCAM_IDPROP
DOUBLE dbVal	Current value needs to be set
INT32 nChn	Current channel needs to be obtained (the default is 0, monochrome camera is 0)

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_IDPROP	Invalid parameter, when product information code does not exist, refer to TUCAM_IDPROP
TUCAMRET_INVALID_VALUE	Invalid value, when pAttr pointer is empty
TUCAMRET_NOT_SUPPORT	When the underlying request is not supported
TUCAMRET_OUT_OF_RANGE	When the channel needs to be obtained exceeds the range
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Prop_GetAttr, TUCAM_Prop_GetValue, TUCAM_Prop_GetValueText

TUCAM_Prop_GetValueText**Description**

Obtain the text information of current attribute value of attribute parameters. The specific supported performance type can be referred to TUCAM_IDPROP.

Statement

```
TUCAMRET TUCAM_Prop_GetValueText (HDTUCAM hTUCam, PTUCAM_VALUE_TEXT pVal);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_VALUE_TEXT pVal	Obtain the text information structure pointer of performance parameters, TUCAM_VALUE_TEXT

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_FAILURE	Text buffer size is 0 or when pText pointer is empty
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Prop_GetAttr, TUCAM_Prop_SetValue, TUCAM_Prop_GetValue

TUCAM_Buf_Alloc**Description**

Allocate memory for data capture. When the application calls this interface, SDK will allocate the necessary internal buffer to buffer image acquisition. Capture does not start from this moment. When the acquisition is started, the application must call [TUCAM_Cap_Start](#) interface. If the buffer is no longer necessary, the application should call [TUCAM_Buf_Release](#) interface to release the internal buffer.

Statement

```
TUCAMRET  TUCAM_Buf_Alloc (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_FRAME pFrame	Picture time frame structure, refer to TUCAM_FRAME

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_PARAM	when pFrame pointer is empty
TUCAMRET_EXCLUDED	When TUCAM_Buf_Alloc is called and not released
TUCAMRET_OUT_OF_RANGE	When the number of frame need to be obtained exceeds the range
TUCAMRET_NO_MEMORY	When memory is low
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame, TUCAM_Buf_CopyFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop

TUCAM_Buf_Release**Description**

Free up memory space for data capture. If the interface is called during capture, the interface will return to the state that the camera is busy.

Statement

```
TUCAMRET  TUCAM_Buf_Release (HDTUCAM hTUCam);
```

Parameters

HDTUCAM hTUCam	Camera handle
----------------	---------------

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_BUSY	Camera busy
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open, TUCAM_Dev_Close
TUCAM_Buf_Alloc
TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame, TUCAM_Buf_CopyFrame
TUCAM_Cap_Start, TUCAM_Cap_Stop

TUCAM_Buf_AbortWait

Description

It is used for the waiting when stopping the data capture. After calling [TUCAM_Buf_WaitForFrame](#) for data capture waiting, use this interface to abort waiting.

Statement

```
TUCAMRET  TUCAM_Buf_AbortWait (HDTUCAM hTUCam);
```

Parameters

HDTUCAM hTUCam	Camera handle
----------------	---------------

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open, TUCAM_Dev_Close
TUCAM_Buf_Alloc, TUCAM_Buf_Release
TUCAM_Buf_WaitForFrame, TUCAM_Buf_CopyFrame
TUCAM_Cap_Start, TUCAM_Cap_Stop

TUCAM_Buf_WaitForFrame

Description

It is used for the completion of data capture. By calling **TUCAM_Buf_Alloc** the space allocated, the captured frame data is obtained. It must be used after calling **TUCAM_Cap_Start** to start capturing; otherwise it will return the state of not ready.

This function belongs to choke function, until the data capture is completed or call **TUCAM_Buf_AbortWait** to abort.

Set the frames need to be captured in uiRsdSize of frame structure, and it is valid for the trigger mode. For example: When 5 frames need to be returned in one trigger, this function will be returned after waiting for the end of capture of five frames of data.

Note: The data ordering of the returned frame structure pBuffer is frame head (usHeader) + image data (uiImgSize) + reserved bits (uiHstSize). If the multi-frame is returned, then it is arranged in this order.

Statement

```
TUCAMRET  TUCAM_Buf_WaitForFrame (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_FRAME pFrame	Frame structure pointer

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NOT_READY	When starting capture without calling TUCAM_Cap_Start
TUCAMRET_NO_MEMORY	When TUCAM_Buf_Alloc is not called to create memory space
TUCAMRET_NO_RESOURCE	When pFrame pointer is empty
TUCAMRET_OUT_OF_RANGE	When the number of frames need to be obtained is greater than 1 and the format obtained is different from TUCAM_Buf_Alloc
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit

TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_CopyFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop

TUCAM_Buf_CopyFrame

Description

It is used for the image format data copied after the completion of waiting data capture different from **TUCAM_Buf_Alloc**. It must be called after **TUCAM_Buf_WaitForFrame** is returned; otherwise the correct image data cannot be obtained.

For example: when the allocated image format is TUFRM_FMT_RGB888, data of other formats can be copied by using this function (i.e., TUFRM_FMT_RAW), and this interface cannot copy the data larger than 1 frame, that is, uiRsdSize in the frame structure cannot be larger than 1.

Note: The data ordering of the returned frame structure pBuffer is frame head (usHeader) + image data (uiImgSize) + reserved bits (uiHstSize). Return of the multi-frame data is not supported.

Statement

Statement

```
TUCAMRET  TUCAM_Buf_CopyFrame (HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_FRAME pFrame	Frame structure pointer

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NOT_READY	when starting capture without calling TUCAM_Cap_Start
TUCAMRET_NO_MEMORY	When TUCAM_Buf_Alloc is not called to create memory space
TUCAMRET_NO_RESOURCE	when pFrame pointer is empty
TUCAMRET_OUT_OF_RANGE	When the number of frames need to be obtained is greater than 1 and the format obtained is different from TUCAM_Buf_Alloc
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop

TUCAM_Cap_SetROI

Description

It is used to set the interested areas of image, with the upper left corner as the origin of coordinates. The set horizontal offset, vertical offset, width and height must be in multiples of 4.

Statement

```
TUCAMRET  TUCAM_Cap_SetROI (HDTUCAM hTUCam, TUCAM_ROI_ATTR roiAttr);
```

Parameters

HDTUCAM hTUCam	Camera handle
TUCAM_ROI_ATTR roiAttr	Object of ROI attribute structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NOT_SUPPORT	ROI settings not supported
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop
 TUCAM_Cap_GetROI

TUCAM_Cap_GetROI

Description

It is used to set the interested areas of image, with the upper left corner as the origin of coordinates. The set horizontal offset, vertical offset, width and height must be in multiples of 4.

Capture does not start from this moment. Acquisition starts after [TUCAM_Cap_Start](#) interface is called.

Statement

```
TUCAMRET  TUCAM_Cap_GetROI (HDTUCAM hTUCam, PTUCAM_ROI_ATTR pRoiAttr);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_ROI_ATTR pRoiAttr	Pointer of ROI attribute structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NOT_SUPPORT	ROI settings not supported
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop
 TUCAM_Cap_SetROI

TUCAM_Cap_SetTrigger

Description

It is used for setting the trigger attribute. Capture does not start from this moment.

Acquisition starts after **TUCAM_Cap_Start** interface is called.

Exposure Mode:

TUCTE_EXPTM means the exposure time is set by the software
 TUCTE_WIDTH means the exposure time is set by the input level width

Stimulation edge mode:

TUCTD_RISING means the trigger signal is rising edge valid
 TUCTD_FALLING means the trigger signal is falling edge valid

Number of frames: it means after receiving a trigger signal, the exposure time of each image is the same no matter how many images are shot, which depends on the software settings. (When choosing level width, the function is invalid.)

Delay: it means after receiving a trigger signal, the delay time of the desired delay time can be set to trigger the camera exposure.




Fig. 19: Normal level trigger mode (rising edge)

Parameters supported by trigger mode:

Mode	TUCAM_TRIGGER_EXP	TUCAM_TRIGGER_EDGE	Delay	Number of frames
Standard trigger	Yes	Yes	Yes	Yes
Synchronization trigger	Yes	Yes	No	No
Global trigger	No	Yes	No	No
Software trigger	No	No	No	No

Synchronization trigger: namely the synchronous picture grabbing, first trigger to start, and the second trigger to output synchronization image.




Fig. 21: Normal synchronous readout trigger mode (rising edge)

Global trigger: generally used for the scenes of controllable light source.




Fig. 20: Global exposure level trigger mode

Software trigger: Simulate trigger signal via the software command.

Statement

```
TUCAMRET TUCAM_Cap_SetTrigger (HDTUCAM hTUCam, TUCAM_TRIGGER_ATTR tgrAttr);
```

Parameters

HDTUCAM hTUCam	Camera handle
TUCAM_TRIGGER_ATTR tgrAttr	Object of trigger attribute structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close

TUCAM_Buf_Alloc, TUCAM_Buf_Release
TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
TUCAM_Cap_Start, TUCAM_Cap_Stop
TUCAM_Cap_GetTrigger, TUCAM_Cap_DoSoftwareTrigger

TUCAM_Cap_GetTrigger

Description

It is used to obtain the trigger attributes.

Statement

```
TUCAMRET  TUCAM_Cap_GetTrigger (HDTUCAM hTUCam, PTUCAM_TRIGGER_ATTR pTgrAttr);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_TRIGGER_ATTR pTgrAttr	Pointer of trigger attribute structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open, TUCAM_Dev_Close
TUCAM_Buf_Alloc, TUCAM_Buf_Release
TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
TUCAM_Cap_Start, TUCAM_Cap_Stop
TUCAM_Cap_SetTrigger

TUCAM_Cap_DoSoftwareTrigger

Description

Execute software trigger commands.

Statement

```
TUCAMRET  TUCAM_Cap_DoSoftwareTrigger(HDTUCAM hTUCam);
```

Parameters

HDTUCAM hTUCam	Camera handle
----------------	---------------

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_FAILURE	Failed to execute trigger commands
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop
 TUCAM_Cap_SetTrigger

TUCAM_Cap_Start**Description**

Start data capture. Prior to capture, the interested areas and trigger mode should be configured.

Statement

```
TUCAMRET  TUCAM_Cap_Start(HDTUCAM hTUCam, UINT32 uiMode);
```

Parameters

HDTUCAM hTUCam	Camera handle
UINT32 uiMode	Camera capture mode

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
-------------------	---------------------------

TUCAMRET_FAILOPEN_BULKIN	Failed to open camera capture
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start
 TUCAM_Cap_SetTrigger, TUCAM_SetROI

TUCAM_Cap_Stop

Description

Stop data capture.

Statement

```
TUCAMRET  TUCAM_Cap_Stop (HDTUCAM hTUCam);
```

Parameters

HDTUCAM hTUCam	Camera handle
----------------	---------------

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_FAILOPEN_BULKIN	Failed to open camera capture
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_AbortWait, TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Stop

TUCAM_File_SaveImage

Description

Save frame data.

Statement

```
TUCAMRET TUCAM_File_SaveImage (HDTUCAM hTUCam, TUCAM_FILE_SAVE fileSave);
```

Parameters

HDTUCAM hTUCam	Camera handle
TUCAM_FILE_SAVE fileSave	File save structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_PARAM	Parameters entered invalid
TUCAMRET_INVALID_PATH	Path entered does not exist
TUCAMRET_FAILURE	Failed to save file
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_WaitForFrame, TUCAM_Buf_CopyFrame

TUCAM_Rec_Start**Description**

Open the video file, and save the video frame data, the data are not written at this time. The frame rate set should be greater than 1fps, and those less than 1fps will be regarded as 1fps to create video files.

Statement

```
TUCAMRET TUCAM_Rec_Start(HDTUCAM hTUCam, TUCAM_REC_SAVE recSave);
```

Parameters

HDTUCAM hTUCam	Camera handle
TUCAM_REC_SAVE recSave	Video file save structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_INVALID_PARAM	Parameters entered invalid
TUCAMRET_INVALID_PATH	Path entered does not exist
TUCAMRET_FAILOPEN_FILE	Failed to open files
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop
 TUCAM_Rec_Stop, TUCAM_Rec_AppendFrame

TUCAM_Rec_AppendFrame**Description**

Write the image data into the file, and call the interface at [TUCAM_Buf_WaitForFrame](#).

Statement

```
TUCAMRET  TUCAM_Rec_AppendFrame(HDTUCAM hTUCam, PTUCAM_FRAME pFrame);
```

Parameters

HDTUCAM hTUCam	Camera handle
PTUCAM_FRAME pFrame	Frame structure pointer

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NOT_READY	TUCAM_Rec_Start interface not called
TUCAMRET_OUT_OF_RANGE	The image width and height are inconsistent with the

ones when created

TUCAMRET_INVALID_CAMERA Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop
 TUCAM_Rec_Start, TUCAM_Rec_Stop

TUCAM_Rec_Stop

Description

Close video file, and calling of [TUCAM_Rec_AppendFrame](#) at this time will not be able to write data.

Statement

```
TUCAMRET  TUCAM_Rec_Stop (HDTUCAM hTUCam);
```

Parameters

HDTUCAM hTUCam Camera handle

Error Codes

TUCAMRET_NOT_INIT TUCAM-API not initialized
 TUCAMRET_INVALID_CAMERA Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
 TUCAM_Dev_Open, TUCAM_Dev_Close
 TUCAM_Buf_Alloc, TUCAM_Buf_Release
 TUCAM_Buf_WaitForFrame
 TUCAM_Cap_Start, TUCAM_Cap_Stop
 TUCAM_Rec_Start, TUCAM_Rec_AppendFrame

TUCAM_Reg_Read

Description

Contents read from the register. Refer to TUREG_TYPE for types of reading.

Statement

```
TUCAMRET TUCAM_Reg_Read (HDTUCAM hTUCam, TUCAM_REG_RW regRW);
```

Parameters

HDTUCAM hTUCam	Camera handle
TUCAM_REG_RW regRW	Register read/write structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NO_MEMORY	The incoming buffer did not allocate memory space
TUCAMRET_NOT_SUPPORT	Reading of this type is not supported
TUCAMRET_INVALID_IDPARAM	Invalid type, refer to TUREG_TYPE
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open, TUCAM_Dev_Close
TUCAM_Reg_Write

TUCAM_Reg_Write

Description

Contents written into the register. Refer to TUREG_TYPE for types of writing.

Statement

```
TUCAMRET TUCAM_Reg_Write(HDTUCAM hTUCam, TUCAM_REG_RW regRW);
```

Parameters

HDTUCAM hTUCam	Camera handle
TUCAM_REG_RW regRW	Register read/write structure

Error Codes

TUCAMRET_NOT_INIT	TUCAM-API not initialized
TUCAMRET_NO_MEMORY	The incoming buffer did not allocate memory space
TUCAMRET_NOT_SUPPORT	Reading of this type is not supported
TUCAMRET_INVALID_IDPARAM	Invalid type, refer to TUREG_TYPE
TUCAMRET_INVALID_CAMERA	Invalid camera, when the camera handle does not exist

Related Interfaces

TUCAM_Api_Init, TUCAM_Api_Uninit
TUCAM_Dev_Open, TUCAM_Dev_Close
TUCAM_Reg_Read